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
- [Modulus™ 微孔板型多功能检测仪](#)
- 吸收光检测模块
- Cortisol EIA Assay Kit (Monobind, Inc., 3625-300), including:
  - **Human Serum references, 1 mL/vial.** Six vials of serum reference for Cortisol at concentrations of 0 (A), 1.0 (B), 4.0 (C), 10.0 (D), 20.0 (E), and 50.0 (F) µg/dL. A preservative has been added. Store at 2 - 8° C.
  - **Cortisol Enzyme Reagent, 1 mL/vial.** One vial of Cortisol (Analog)- horseradish peroxides (HRP) conjugate in a protein-stabilizing matrix with blue dye. Store at 2 - 8° C.
  - **Cortisol Conjugate Buffer, 7 mL/vial.** One vial of reagent containing buffer, red dye, preservative, and binding protein inhibitors. Store at 2 - 8° C.
  - **Cortisol Biotin Reagent, 7 mL bottle.** One bottle of reagent containing anticortisol biotinylated mIgG conjugate in buffer, green dye, and preservative. Store at 2 - 8° C.
  - **Streptavidin-Coated Plate, 96-well plate.** One 96-well microplate coated with 1.0 µg/mL streptavidin and packaged in an aluminum bag with drying agent. Store at 2 - 8° C.
  - **Wash Solution, 20 mL vial.** One vial containing a surfactant in buffered saline. A preservative has been added. Store at 2 - 8° C.
  - **Substrate A, 7 mL/vial.** One vial containing tetramethylbenzidine (TMB) in buffer. Store at 2 - 8° C.
  - **Substrate B, 7 mL/vial.** One vial containing hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) in buffer. Store at 2 - 8° C.
  - **Stop Solution, 8 mL/vial.** One vial containing a strong acid (1N HCL). Store at 2 - 30° C.

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Streptavidin<sub>CW</sub> = Streptavidin immobilized on well

Immobilized complex = sandwich complex bound to the solid surface

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## 6. REFERENCES

1. Burtis, C.A., Ashweed, E.R. Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company, 2nd edition, 1994: 1825 - 1827.
2. Foster, L., Dunn, R. "Single antibody technique for radioimmunoassay of cortisol in unextracted serum or plasma." Clin. Chem.20 (1974): 365.
3. Williams Textbook of Endocrinology. Edited by Wilson, J.D., Foster, D.W. Philadelphia: WB Saunders, 7th edition, 1985.
4. Monobind Data Sheet 3625-300: Monobind Cortisol EIA Kit.